

Medical Refrigerator(2-8°C)



Medical Refrigerator(2-8°C)



Medical Refrigerator(2-8°C)

Specification

Model		MR-5V250	MR-5V310	MR-5V315	MR-5V320	MR-5V415
Capacity (L)		226	312	316	316	416
Size (WxDxH)	Internal (mm)	534x406x1150	560x485x1350	540x425x1380	540x425x1380	685x429x1380
	External (mm)	625x594x1707	650x624x1902	640x554x1925	660x555x1915	785x552x1927
Temperature	Range (°C) at RT. 10 ~ 32°C	2 ~ 8	2 ~ 8	2 ~ 8	2 ~ 8	2 ~ 8
	Sensor	NTC	NTC	NTC	NTC	NTC
	Controller	Microprocessor	Microprocessor	Microprocessor	Microprocessor	Microprocessor
Cooling System	Refrigerant	R600a (CFC Free)	R290 (CFC Free)	R290 (CFC Free)	R290 (CFC Free)	R290 (CFC Free)
	Defrost	Auto	Auto	Auto	Auto	Auto
	Compressor(Qty)	1	1	1	1	1
Display		Digital display	Digital display	Digital display	Digital display	Digital display
Audible &Visible Alarm	High / Low temperature	Standard	Standard	Standard	Standard	Standard
	High ambient temperature	/	/	Standard	/	Standard
	Abnormal voltage	/	/	/	/	/
	Power failure backup system (alarm)	8h	8h	72h	8h	72h
	Door ajar	Standard	Standard	Standard	Standard	Standard
Material	Internal	HIPS	HIPS	Electro-galvanized steel with antimicro-bial powder coating	Electro-galvanized steel with antimicrobial powder coating	Electro-galvanized steel with antimicrobial powder coating
	External	Electro-galvanized steel with antimicro-bial powder coating	Electro-galvanized steel with antimicrobial powder coating	Electro-galvanized steel with antimicro-bial powder coating	Electro-galvanized steel with antimicrobial powder coating	Electro-galvanized steel with antimicrobial powder coating
Accessories	USB port	Optional	/	Standard	Optional	Standard
	Casters	Standard	Standard	Standard	Standard	Standard
	Chart recorder	/	/	/	/	/
	Data logger	Optional	Optional	Optional	Optional	Optional
	Remote alarm port	Standard	Standard	Standard	Standard	Standard
	Test hole (Qty) / Diameter (mm)	1 / 25	1 / 25	1 / 25	1 / 25	1 / 25
Shelves / Drawers (Qty)		Shelves / 4	Shelves / 5	Shelves / 6, Basket / 1	Shelves / 5	Shelves / 6 , Basket / 1
Sound Emission (dB)		50	50	50	50	50
Electricity	Voltage / Frequency	AC110 / 220V±10%, 50 / 60Hz	AC110 / 220V±10%, 50 / 60Hz	AC110 / 220V±10%, 50 / 60Hz	AC110 / 220V±10%, 50 / 60Hz	AC110 / 220V±10%, 50 / 60Hz
	Consumption (W)	140	245	294	248	330
N.W. / G.W. (kg)		73 / 76	86 / 90	92.5 / 100	92.5 / 100	120 / 133
Shipping Size (WxDxH) (mm)		707x676x1745	710x680x1960	710x680x1960	710x680x1960	854x646x2100