

XC-400 Simulator Series



XC-407 Human Trachea Intubation Model

This model is designed as a teaching model for training the medical students, basic clinic nurses and first aid people to demonstrate and practice trachea intubation through mouth.

- 1) Trachea Intubation through mouth.
- 2) Compare pupils: one is normal, another is large.

Packing: 1pcs/carton, 58x37x29cm, 7kgs

XC-407A Advanced Human Trachea Intubation Model

Features:

- 1) Upgraded from XC-407.
- 2) Neck backward function, more realistic simulation of endotracheal intubation technique and teaching demonstration and practice operation.

Packing: 1pcs/carton, 58x37x29cm, 7kgs



XC-410 The Model of Newborn Intubation

This model is used to demonstrate the newborn baby intubation and train the students to learn intubation skills. Served as an ideal model for medical colleges, nursing schools and hospitals etc.

Packing: 12pcs/carton, 54x40x50cm, 10kgs

XC-410A Infant Intubation Training Model

Realistic anatomical structure of infant mouth, pharynx, trachea. Tracheal intubation, directly observe the expansion of lungs and stomach by insufflating air into the tube and stomach so as to verify if the tube position is right.

Packing: 1pcs/carton, 38x20x28cm, 3kgs

