

XC-400 Simulator Series

XC-404 Basic CPR Training Model (Half Body)

The model is designed according to human anatomical half body, served as a teaching aid for education, hygiene, first aid training and demonstration. It shows the half body from the head to the sixth rib of the chest. The cross-section is color painted to show different organs clearly and vivid. The full shape of the heart and lung showing the various movements when demonstrating heart pressure and mouth to mouth respiration.

Feature:

CPR manikin features:

- 1) Head can be raised up and down, neck joint is movable.
- 2) Eyelid can open and close.
- 3) Nose wing is soft, easy to tweak. The mouth can open and close. Can do mouth to mouth respiration or mouth to nose respiration.
- 4) Consist of spring and heart model in heart part.
- 5) Skin is springy with natural skin color.

CPR Instrument feature:

- 1) Green light bright, heart pressure correct.
- 2) Red light bright, heart pressure too strong.
- 3) No light bright, pressure on a wrong place or pressure is too light.

Packing: 1pcs/carton, 65x27x47cm, 8kgs



XC-404A/404B Electronic Half-body CPR Training Manikin (Male/Female)

This Adult CPR model simple, economical and portable, it is used at all levels of hospitals, medical schools, an ideal model of CPR training and teaching.

Feature:

- 1) Airway opening simulation.
- 2) Handling press on the chest: strength correct (4-5cm), show correct bumping. Pressure too large (over 5cm), show alert bumping.
- 3) Mouth to mouth (blowing): the tidal volume range judged by observing the ups and downs of the chest. (tidal volume standard $\leq 500\text{ml}/600\text{ml}-1000\text{ml}$).
- 4) Operating frequency: the latest international standard: 100 times/min.

Packing: 1pcs/carton, 78x36x25cm, 7kgs