



Wincom Company Ltd.

Wind Blow Digital Thermostat Laboratory Oven

Model NO. DHG-series

Certification: ISO,

Structure: Floortop,

Material : Stainless Steel,

Type : Oven,

Transport Package: Foam+Carton,

Specification : 1-9999min.

Trademark : OEM, Origin

Brand : Wincome

Made in : China



Product Description

Product Description:

1. Approval quality.
2. In Neutro-packing
3. Certification: ISO.

Wind Blow THERMOSTAT OVEN		
Model No	Working room size(mm)	Remark
DHG-9023A	300x300x280(25L)	• Timer range: 1-9999min. • Temp.range: RT+5-250°C • Temp.fluctuation :±1°C • A: Stainless steel Workroom • S: The temperature controller with digital display • With Fan for Wind blow
DHG-9053A	415x370x345(53L)	
DHG-9123A	554x370x524(100L)	
DHG-9101-0SA	350X350X350(42L)	
DHG-9101-1SA	350X450X450(70L)	
DHG-9101-230SA	550X450X550(225L)	