

Sieve Shaker Device

Model : RTC-200/300

Manufacture ring Country: Bangladesh.

Test sieve shaker for the fields of pharmacy, metallurgy, foods, cosmetics and grain, and characterized by convenient operation, no rotary component, and electronic controlled amplitude and frequency.

Material	Stainless steel
Layers	1-9 Layers
Diameter	0.25-3mm
Noise	≤50dB
Input quantity	≤200
Amplitude	≤5mm
Voltage	220V 50Hz
Speed	1450r/m
Power	0.125 kw
Sieves frame	430 x 430 x 750 mm
Weight	40Kg
Sieves USE	200mm, 300mm
Motor	Gazi, Bangladesh
Warranty	3 Years

