

## XC-400 Simulator Series



### XC-439 Advanced Suture Practice Arm



### XC-440 Advanced Suture Practice Leg

These two products provide a realistic suturing demonstration. Soft and pliable for easy sewing. These products are provided with "wounds" which can be sutured repeatedly, new "wounds" can be cut also.

XC-439 Packing: 4pcs/carton, 64x33x29cm, 8kgs

XC-440 Packing: 2pcs/carton, 74x43x29cm, 10kgs

### XC-471 Surgical Suture and Bandaging Model

The simulator provides commonly used surgical incisions, which simulate all kinds of female thoracic and abdominal incisions, tube drains and amputation stump. Students can practice wounds treatments of washing, disinfecting, changing dressing, bandaging and etc.

Simulator features the following simulated wounds:

1. Thyroidectomy
2. Mid-sternal split - with two chest tube drains
3. Right mastectomy - with simulated drain
4. Right cholecystectomy - with simulated T-tube
5. Laparotomy.
6. Appendectomy.
7. Left colostomy.
8. Right ileostomy.
9. Abdominal Hysterectomy .  
(transverse incision).
10. Left nephrectomy .  
(thoraco-abdominal incision).
11. Right nephrectomy (oblique incision).
12. Laminectomy.
13. Sacral decubitus ulcer - stage 2.
14. Right thigh amputation stump.

Packing: 1pcs/carton, 79x48x27cm, 9kgs

